

Spare Parts List Evita 4 / 2 dura, US version




Table of contents:


Picture	Description
	PM KITS, SENSORS, ACCESSORIES, SOFTWARE KITS
2, 3	EVITA 2 DURA
4, 5	EVITA 4
6 - 11	PNEUMATICS, EVITA 4 / 2 DURA
14	TUBING PLAN
15	BLOCK BETWEEN MIXER AND PNEUMATIC BASE PLATE
16	PLATE BETWEEN MIXER AND PNEUMATIC BASE PLATE
17	HEATING SUPPORT FOR EXP. VALVE
18	INSPIRATION BLOCK, CPL.
19	PARALLEL MIXER
20	GAS SUPPLY BLOCK (DRAEGER + FAS)
21	PRESSURE METERING BLOCK, CPL.
22	PATIENT SYSTEM (EXP. VALVE)
12	OPTION SPO2
13	OPTION CO2
23, 24	ELECTRONICS (PCBs) EVITA 4 / 2 DURA
25	FRONT PANEL ELECTRONICS, EVITA 4
26 - 28	OPERATING DEVICE, EVITA 4
29	DISPLAY TOSHIBA, VERSION 1, EVITA 4
30	DISPLAY TOSHIBA, VERSION 0, EVITA 4
31	DISPLAY SHARP, VERSION 3, EVITA 4
32	DISPLAY NEC, VERION , EVITA 4
65, 66	FRONTPANEL, EVITA 2 DURA
67	MODIFICATION KIT LC-DISPLAY, EVITA 2 DURA
68	MODIFICATION KIT COLOR TFT-DISPLAY, EVITA 2 DURA
33	HINGED SUPPORT ARM
34 - 44	TROLLEY, VERSION WITH 2 COLUMNS
34	SUPPORT
35, 36	CARRIAGE AND BOTTOM SECTION

Picture	Description
37, 38	RECEIVER WITH SUPPORT AND FLAP
40	SIDE PART, RIGHT
41	SIDE PART, LEFT
44	FIXING SET FOR CUPBOARD EVITA
59 - 64	TROLLEY, VERSION WITH 1 COLUMN, CURRENT VERSION
59	TROLLEY WITH 1 COLUMN
60	TROLLEY WITH 1 COLUMN
61	OPTION HUMIDIFIER RACK.
62	OPTION CYLINDER SUPPORT
64	KIT BATTERY CABLE EVITA MOBIL

PM KITS, SENSORS, ACCESSORIES, SOFTWARE KITS


Description	Order-Code
PMS Kits For more information (parts needed and replacement frequency) refer to the PMS Instructions Evita 4 / 2 dura.	
Kit PMS 12 months - Evita 4 / 2dura with 2 Expiration Valves	4114761
2 x Filter	8412384
2 x Diaphragm	8410181
2 x Sealing washer	8407979
Kit PMS 24 months - Evita 4 / 2dura with 2 Expiration Valves	4114762
1 x Kit PMS 12 month	4114761
2 x Filter	8408208
1 x Flat packing gasket (for Dräger gas connection block)	8408204
2 x O-Ring (for FAS gas connection block)	8411516
1 x Lithium Battery	1835343
Note: There are no batteries for the DC option included in this kit, refer to the PMS Instructions Evita 4 / 2 dura.	
SENSORS AND SENSOR ACCESSORIES, for more information refer to the Operating Instructions	
O2 Sensor	6850645
Spirolog sensor (exp. flow measurement), 5 pack	8403735
CO2 sensor	6870300
Adult cuvette	6870279
Pediatric cuvette	6870280
Neonatal flow sensor ISO 15 (sensor housing with 1 flow sensor element)	8411130
Neonatal flow sensor element, 5 pack	8410179
Connecting cable neonatal flow sensor	8409626
ACCESSORIES, for more information refer to the Operating Instructions	
Evita AIR filter kit	4510808
AIR supply filter and water trap, straight	D700133
Evita O2 filter kit	4510807
O2 supply filter and water trap, straight	D700131

Description	Order-Code
SOFTWARE KITS, for more information refer to the Operating Instructions	
Evita 4, SW 4.00 (update)	8414467
Evita 2 DURA, SW 4.00 (update)	8414468
TEST LUNGS	
Test lung (blue, Draeger)	8403201
Test lung (white)	8401892
Test lung, pediatric	8409742

EVITA 2 dura**Picture 2**

EVITA 2 dura


Picture 3




Item No.	Description	Order-Code
1	VENT PLUG	8412857
2	VENT PLUG	8412858
3	SCREW	2M17288
4	POWER CORD US	1841793
5	SCREW	8412378
6	EXPIRATION VALVE, SEE PICTURE 22	
7	KNURLED SCREW	8306225
8	VENT PLUG	G60455
9	BAR	8406763
10	CAP, WHITE	M25146
11	FILTER	8412384
12	PARKING SUPPORT	8411784
	HOOD 4H, 1B, WITHOUT ILLUSTRATION NOTE: ORDER IN ADDITION LABEL P/N 8306124	8412379
	LABEL "DANGER, WARNING, CAUTION" FOR HOOD	8306124

EVITA 4

Picture 4





EVITA 4**Picture 5**

Item No.	Description	Order-Code
2	VENT PLUG	8412857
3	VENT PLUG	8412858
4	HOOD 4H, 1B NOTE: ORDER IN ADDITION LABEL P/N 8306124	8412379
4A	LABEL "DANGER, WARNING, CAUTION" FOR HOOD	8306124
5	SCREW	2M17288
6	OPERATING DEVICE (EVITA 4) SEE PICTURE 26 - 32	
7	FRONT PANEL, CPL., SEE PICTURE 25	
8	ELECTRONIC, SEE PICTURE 23 - 25	
10	SUPPLY MAIN US	1841793
11	CABLE, EVITA 4 - OPERATING DEVICE	1839780
12	SCREW	8412378
13	HOOD	8412375
14	EXPIRATION VALVE , SEE PICTURE 22	8410580
15	PNEUMATIC, SEE PICTURE 6 - 11	
16	KNURLED SCREW	8306225
17	SIDE PANEL 2 H	8408369
18	SIDE WALL 2H, RIGHT	8412355
19	SOCKET	8406504
20	VENT PLUG	G60455
21	BAR	8406763

PNEUMATICS (GAS SUPPLY BLOCK, HOUSING), EVITA 4 / 2 DURA


Picture 6



PNEUMATICS (PCBs, BOARD GUIDE), EVITA 4 / 2 DURA**Picture 7**


**PNEUMATICS (PEEP VALVE, PRESSURE MEASUREMENT BLOCK),
EVITA 4 / 2 DURA**


Picture 8




**PNEUMATICS (INSPIRATION BLOCK, PARALLEL MIXER),
EVITA 4 / 2 DURA**

Picture 9



PNEUMATICS (SIDE PART 2H), EVITA 4 / 2 DURA**Picture 10**


PNEUMATICS (SIDE PART 2H), EVITA 4 / 2 DURA**Picture 11**

Item No.	Description	Order-Code
1-53	PNEUMATIC	8412360
1	BASE PLATE	8412376
2	RUBBER BUFFER	2M04819
3-4	COVER PLATE BACK PANEL, CPL.	8412718
3	DISSIPATOR	8412383
5	CABLE CLIP (2,4X92)	8712007
6	SOCKET	8412361
7	SIDE PART 2H	8412354
8	CABLE BUSHING	8412377
9	CHASSIS PNEUMATIC	8412339
10	KNURLED SCREW	8306225
11	CABLE HARNESS HPSV	8306545
12	PCB PNEUMATIC-MOTHERBOARD	8306571
13	GAS SUPPLY BLOCK, SEE PICTURE 20	
14	BOARD GUIDE	8600527
15	PCB PNEUMATIC CONTROLLER	8306601
16	PCB HPSV-CONTROLLER	8306561
17	SUPPORT, CPL.	8412373
18	SOCKET	8408312
19	PLATE	8412341
20 - 21	PEEP/PIP VALVE COMPLETE WITH CABLE	4116593
20	PEEP/PIP VALVE	8410717
21	CABLE HARNESS PEEP-VALVE	8306557
22	ANGLE BRACKET, CPL.	8412369
23	PRESSURE METERING BLOCK, SEE PICTURE 21	
24	PARALLEL MIXER, SEE PICTURE 19	8411550
25	LABEL	8412803
26	LABEL	8412802
28	INSULATING PLATE	8412211
30	SOCKET	8407712

Item No.	Description	Order-Code
31	PLATE	8412605
32	HEATING SUPPORT, SEE PICTURE 17	
33	CABLE HARNESS FLOWSENSOR	8306549
34	HOSE 4X1,5-SI 50 SH A NF	1190520
35	CONNECTOR RECEIVER	8410663
36	DISC	8300841
37	CABLE HARNESS FLOW-SWITCH	8306556
38	RUBBER PLATE	8410917
39	THREADED SOCKET	8408348
40	SUPPORT	8412365
41	SOCKET	8412362
42	FAN PNEUMATIC, CPL.	8306560
43	PLATE (EVITA 4)	8412222
44	LIP SEAL	8407689
45	CONNECTING HOUSING	8410570
46	SPRING	8410631
47	INSPIRATION BLOCK, SEE PICTURE 18	
48	HINGE, CPL.	8412707
49	SHEET METAL, CPL.	8412382
50	CATCH	2M17303
51	LINK	8412363
52	PIN	8408902
53	FIXING STRAP	8406526
	BASE, without illustration	2M17302
	BUFFER, without illustration	6800478
	HEADER PROFILE 1-2CH15956 SW, without illustration	1191640
	ADH.TAPE 10X4-PUE-FOAM SW, without illustration	1209310

OPTION SPO2, EVITA 4 / 2 DURA


Picture 12



Item No.	Description	Order-Code
1-10	KIT EVITA 4 SAT	8413035
1	DURASENSOR DS-100A	8201001
2	EXTENSION CABLE SENSOR 2,6M	8600859
3	SCREW M3X5 PA	8411753
4	PACKING RING	D04766
5	CABLE HARNESS SPO2	8600787
6	CABLE CLIP (2,4X92)	8712007
7	SAO2 MODULE	8600481
8	CABLE HARNESS SPO2 POWER	8411728
9	DC/DC TRANSFORM. +-15V LV01895	1837869
10	FLAP/MOUNTING SPO2-CABLE	8600677

OPTION CO2, EVITA 4 / 2 DURA


Picture 13



Item No.	Description	Order-Code
	OPTION CAPNO PLUS	8411720
1	PCB CO2-CARRIER, SEE PICTURE 23	
2	CABLE HARNESS CO2-POWER RS 232	8306567
3	PCB POWER	8350411
4	PCB SIGNAL PROCESSOR (CO2-MAIN)	6870296
4a	SOFTWARE 1.13 (CO2-MAIN)	6870508
5	PACKING RING	D04766
6	DISTANCE BOLT	8305182
7	CABLE HARNESS CO2 SENSOR	8306766
	CO2 CUVETTE ADULT, without illustration	6870279
	CO2 CUVETTE PAED, without illustration	6870280
	CO2 MAINSTREAM SENSOR, without illustration, Dräger	6870300
	REPAIR EXCHANGE MAINSTREAM SENSOR (6870300), without illustration, Dräger	6870342
	TEST FILTER, without illustration	6870281
	PARKING SUPPORT, CPL.	8412840
	XFL SHORT CIRCUIT SOCKET 2POL., without illustration	1814508
	SET OF INSCRIPTED CLIPS, without illustration	6870343

TUBING PLAN, EVITA 4 / 2 DURA


Picture 14




Item No.	Description	Order-Code
1	HOSE 4X1 PAE WS	1210165
2	HOSE 4X1PAE SW/WS M31989	1210203
3	HOSE 2X1,5 SI NF 8403323	1203622
4	HOSE 2X1-SIGN GREEN	1204807
5	HOSE 2X1-SIGN ORANGE	1204823
6	HOSE 2X1-SIGN RED	1204793
7	HOSE 2X1-SIGN BROWN	1204815
8	HOSE 2X1-SIGN BL	1204785
9	HOSE 4X1,5-SI 50 SH A NF	1190520
10	TUBE	2M17645
11	T-SOCKET	6800187
12	T-PIECE	8401083
13	DOSAGE	8408222
14	DOSAGE	8412201
15	DOSAGE	8412998

**BLOCK BETWEEN MIXER AND PNEUMATIC BASE PLATE,
EVITA 4 / 2 DURA**


Picture 15




Item No.	Description	Order-Code
1	SHEET METAL	8409177
2	MICROSWITCH	8407265
3	HOUSING (HPSV CONNECTION)	8412976
4	O-RING SEAL	2M08777
5	O-RING SEAL	E20567
6	SINTERED FILTER	8407827
7	SOCKET	8408197

**PLATE BETWEEN MIXER AND PNEUMATIC BASE PLATE,
EVITA 4 / 2 DURA****Picture 16**

Item No.	Description	Order-Code
1	O-RING SEAL	2M08777
2	DISTANCE PLATE	8412212
3	SINTERED FILTER	8407827

HEATING SUPPORT FOR EXP. VALVE, EVITA 4 / 2 DURA**Picture 17**

Item No.	Description	Order-Code
	HEATING SUPPORT, COMPLETE WITH HEATING RESISTORS AND FAN	8412734
4	XLPSK BUL2 AWG24 LV01882	1838318
5	CABLE CLIP (2,4X92)	8712007
	FAN 12 VDC, 40X40X10, without illustration Note: used in units shipped prior to 2000, check voltage labeling on fan before ordering, order fan with the same voltage labeling, P/N 1837702 (12V) or P/N 1844709 (24V)	1837702
	FAN 24 VDC, 40X40X10, without illustration Note: used in units shipped as of 2000, check voltage labeling on fan before ordering, order fan with the same voltage labeling, P/N 1837702 (12V) or P/N 1844709 (24V)	1844709


INSPIRATION BLOCK, CPL. , EVITA 4 / 2 DURA
Picture 18


Item No.	Description	Order-Code
1-32	INSPIRATION BLOCK, COMPL. NOTE: ORDER IN ADDITION LABEL P/N 8306123	8412981
1-22, 28-32	INSPIRATION BLOCK (W/O. P/N 8411749, COVER, CPL. = O2 AMPLIFIER) NOTE: ORDER IN ADDITION LABEL P/N 8306123	8411843
	LABEL "Do not block air intake", FOR INSPIRATORY BLOCK	8306123
1	COVER 1	8412955
2	WASHER	8400812
3	DIAPHRAGM	8403352
4	RECEIVER	8404257
5	CLOSING SPRING	D13509
6	O-RING SEAL	R27696
7	RECEIVER	8412951
8	O-RING SEAL	E20567
9	O-RING SEAL	2M08777
10	PRESSURE CONTROL VALVE	8411745
11	O-RING SEAL	R26807
13	SEALING WASHER	8410307
14	SPRING	2M12034
15	SCREW	8412952
16	CHOKE SCREW	E08481
17	O-RING SEAL	R31296
18	CLOSING BOLT	8412953
19	COVER 2	8412954
20	LIP SEAL	8412938
21	PCB O2-CONTACT	8306671
22	GASKET	8413031
23-27	COVER, CPL. (O2 AMPLIFIER)	8411749
24	DIAPHRAGM	8410682
28	LUER-LOCK, FEMALE	8600213

Item No.	Description	Order-Code
29	SOCKET	8412939
30	CONNECTION PIECE	8410557
31	O-RING SEAL	M20622
32	ANGULAR PORCELAIN BUSH	8410676

PARALLEL MIXER, EVITA 4 / 2 DURA


Picture 19




Item No.	Description	Order-Code
1-14	PARALLEL MIXER	8411550
2,4,8-10	VALVE AIR	8412128
3,4,8-10	VALVE O2	8412126
4	ADHESIVE LABEL	8302721
5	KNURLED SCREW	8306225
6	REP. SET HIGH PRESSURE SENSOR	8306248
7	O-RING SEAL	8410713
8	TOROIDAL SEALING RING	R16442
9	O-RING SEAL	R22344
10	O-RING SEAL	E09548
11	THRUST COLLAR 6, WHITE	M31603
12	THRUST COLLAR 6, BLACK	M31601
13	ANGLE CONNECTION	M30961

GAS SUPPLY BLOCK (DRAEGER + FAS), EVITA 4 / 2 DURA

Picture 20




Item No.	Description	Order-Code
1-18	GAS SUPPLY BLOCK (Dräger)	
3	VALVE CONE	D21138
5	FLAT PACKING	8408204
6	FILTER (FOR P/N 8411848)	8408208
7	RECEIVER	8412871
8	PRESSURE REGULATOR	8408205
9	SOCKET	8408197
10	MICRO-ELECTROVALVE	8412209
11	PLUG-TYPE CONNECTION	M30960
12	THRUST COLLAR 6, WHITE	M31603
13	THRUST COLLAR 6, BLACK	M31601
14	ADHESIVE LABEL	8412928
15	CONNECTION AIR (DISS)	8412962
17	DISK	M19311
18	ADAPTER O2 (DISS)	8412961
	GAS SUPPLY BLOCK, without illustration (FAS)	8411848
	DIAPHRAGM FPM (FOR 8411848), without illustration	8411513
	VALVE TAPPET (FOR 8411848), without illustration	8411514
	NONRETURN VALVE (FOR 8411848), without illustration	8411515
	O-RING SEAL (FOR 8411848), without illustration	8411516
	FILTER (FOR 8411848), without illustration	8408208
10	MICRO-ELECTROVALVE	8412209
15	ADAPTER AIR (DISS)	8412962
17	DISK	M19311
18	ADAPTER O2 (DISS)	8412961

PRESSURE MEASUREMENT BLOCK, EVITA 4 / 2 DURA**Picture 21**

Item No.	Description	Order-Code
1-10	PRESSURE METERING BLOCK, CPL.	8412210
2	ELECTROVALVE	8412993
4	PACKING RING	D04766
5	DISC	8305747
6	PAW SENSOR	8305623
7	INSULATING BUSH	6804141
8	O-RING SEAL	8410713
9	SOCKET	8408197
10	CABLE HARNESS PRESSURE SENSOR	8306538

PATIENT SYSTEM (EXP. VALVE) , EVITA 4 / 2 DURA


Picture 22



Item No.	Description	Order-Code
1-12	EXPIRATION VALVE	8410580
1	VALVE BOX	8410571
2-3	DIAPHRAGM	8412015
2	SEALING WASHER	8407979
3	DIAPHRAGM	8410181
4	O-RING SEAL	E22651
5	VALVE COVER	8410572
6	DIAPHRAGM	8410682
7	GASKET	8410576
8	TUBE	8410574
9	CONNECTING CAP	8410672
10	O-RING SEAL	M29902
11	ANGULAR PORCELAIN BUSH	8410575
12	SIEVE	8412202
13	MODIFICATION KIT WATER TRAP	8413125
13a	POT	8403976


ELECTRONICS (PCBs) , EVITA 4 / 2 DURA

Picture 23



ELECTRONICS (SIDE PART 4H) , EVITA 4 / 2 DURA

Picture 24




Item No.	Description	Order-Code
1	PCB MOTHERBOARD	8306581
2	KNURLED SCREW	8306225
3	CABLE HARNESS PNEUMATIC	8306546
5	ANGLE, CPL.	8412367
6	EMI-SEALING	8600784
7	BOARD GUIDE	8600527
8	LINK	2M17310
9	LINK	8412363
10a	AC/DC POWER PACK COMPLETE	4115870
10b	AC POWER PACK	8306520
10c	RETROFITTING KIT DC-MODULE	8413034
10d	DC POWER PACK (SPARE PART ONLY)	8306530
10e	BATTERY GSBL 12V/3,5AH, ORDER 2 PIECES (formerly p/n 1845284)	1851926
10d	FAN FOR AC POWER PACK	8411886
10e	POWER PACK FUSE (not pictured) (available August 2003)	1850377
11	SIDE PART 4H	8304117
12	DUMMY PLATE	8306547
13	FILTER HOOD	8412385
14	FILTER	8412384
15	PCB CPU 68332	8306591
16	LITHIUM STOR.BATTERY 3V/1400	1835343
16a	REAL-TIME-CLOCK	1837087
17A	PCB CO2-CARRIER WITH BUZZER ON BORAD, SEE TSB EVITA 4 # 39 OR EVITA 2 DURA # 19 FOR FURTHER INFORMATION	8306611
17B	KIT CO2 CARRIER PCB WITH SEPARATE BUZZER, AS OF SW 3.22, SEE EVITA 4 # 39 OR EVITA 2 DURA # 19 FOR FURTHER INFORMATION	8414362
18	CHASSIS ELECTRONIK	8412374
19	TUB, CPL.	8412706

Item No.	Description	Order-Code
20	CABLE HARNESS	8306559
21	XDSUB LOCK BOLT	1822357
22	HEADER PROFILE 1-2CH15956 SW	1191640
23	SUPPORT	2M17304
25	OPTION CAPNO PLUS, SEE PICTURE 13	8411720
	SPO2 MODULE, without illustration	8600481
	PCB PAEDIATRIC-FLOW, without illustration	8306631
	PCB COMMUNICATION, without illustration	8306621

FRONT PANEL ELECTRONICS, EVITA 4


Picture 25



Item No.	Description	Order-Code
1	FRONT PANEL	8412639
2	BINDER	8412884
3	XDSUB LOCK BOLT LV00575 consisting of: bolt, nut, serrated lock washer	1822357
4	CABLE HARNESS,FRONT (EVITA 4)	8306775
5	FRONT PLATE, CPL.	8412839


OPERATING DEVICE, EVITA 4

Picture 26




OPERATING DEVICE (LOCK) , EVITA 4

Picture 27



OPERATING DEVICE (SHAFT ENCODER) , EVITA 4

Picture 28




Item No.	Description	Order-Code
1-41	OPERATING DEVICE (EVITA 4)	8412720
1-41	REP. EXCH. OPERAT. DEVICE (8412720)	8411683
1	RAIL	8412627
2	UNLOCK KEY	8412632
3	SPRING	M14362
4	GUIDING PIN	8412625
5	SPRING	M09360
6	INSERT NUT	8412626
7	HOUSING (HINGE)	8412611
8	CLOSING SPRING	M14164
9	LOCK	8412622
11	ROLLER	8412631
12	BUFFER	6800478
13	XDSUB LOCK BOLT	1822357
14	CABLE HARNESS TOUCH-SCREEN	8412898
15	HOUSING (OPERATING DEVICE)	8412601
16	SPACER	8412808
17	SPRING STEEL SHEET	8412703
18	SQUARE NUT	8412883
19	CABLE HARNESS, KEY PAD	8412897
20	PLATE	8412885
21	PCB GRAPHIC CONTROLLER	8306641
22	PACKING RING	D04766
23	SHEET METAL, CPL.	8412843
24	TOGGLE	8412637
25	SPRING	E24916
26	SOCKET	8412844
27	BINDER	8412884
28	REPAIR KIT DISPLAY (EVITA 4) INCLUDING ALL PARTS AND THE LATEST DISPLAY VERSION	8411711

Item No.	Description	Order-Code
	LCD BACK LIGHTING:	
	TOSHIBA, version 1 ('95), SEE PICTURE 29	
	NEC, version 2 ('95, '96), SEE PICTURE 32	
	SHARP 3, version 3 ('96, '97), SEE PICTURE 31	
	TOSHIBA, version 0, as of '98, SEE PICTURE 30	
29-30	LOUDSPEAKER, CPL.	8306558
30	XLPSK BUL2 AWG24	1838318
31	TOUCHSCREEN CPL. WITH PCB	8306638
31a	FRAMED SCREEN, CPL. (EVITA 4), W/O PCB	8411746
33	O-RING SEAL	8412941
34	DISPLAY, CPL.	8412137
35	OUTER FRAME	8412602
36A	CONTROL KNOB	M29655
36B	LEAF SPRING (FOR CONTROL KNOB)	M29952
37	BUSH	8412893
38	SHAFT ENCODER, CPL.	8306565
39	ADHESIVE LABEL	8303296
40	PLUG-IN STRIPE (INTERN.)	8412380
41	SILICON GASKET 1,5X0,5 DIN 7715-EZ SI/NF	1212184

DISPLAY TOHIBA, VERSION 1, EVITA 4


Picture 29




Item No.	Description	Order-Code
	REPAIR KIT DISPLAY (EVITA 4) INCLUDING ALL PARTS AND THE LATEST DISPLAY VERSION	8411711
1	MULTIPLEXER COLOUR TOSHIBA	8306640
2	LOOKUP TFT-TOSHIBA	8306795
5	PACKING RING	D04766
	LCD BACK LIGHTING(TOSHIBA 1) VERSION 1, without illustration	1841424

DISPLAY TOHIBA, VERSION 0, EVITA 4

Picture 30




Item No.	Description	Order-Code
	REPAIR KIT DISPLAY (EVITA 4) INCLUDING ALL PARTS AND THE LATEST DISPLAY VERSION	8411711
1	LC-DISPLAY 640X480	1843540
2	PIXELMUXER TOSHIBA LTM 10C209A	8411859
3	LOOKUP TFT-TOSHIBA	8306795
4	U DC/AC-CONVER 12V/550	1843559
5	PACKING RING	D04766
6	SUPPORT, RIGHT	8411862
7	SUPPORT, LEFT	8411861
8	GASKET	D20682
9	CABLE CLIP	8712065
	LCD-BACK LIGHTING TOSHIBA EVITA4, VERSION 0, WITHOUT ILLUSTRATION	1844776

DISPLAY SHARP, VERSION 3, EVITA 4**Picture 31**

Item No.	Description	Order-Code
	REPAIR KIT DISPLAY (EVITA 4) INCLUDING ALL PARTS AND THE LATEST DISPLAY VERSION	8411711
1	MULTIPLEXER SHARP	8411776
2	LOOKUP TFT-TOSHIBA	8306795
8	INVERTER SHARP	8411732
10	SPACER	8300428
11	CABLE CLIP	8712065
12	SCREW M3X5 PA	8411753
	BACK LIGHTING (SHARP), without illustration, version 3 (1996/1997)	1843095

DISPLAY NEC, VERSION 2, EVITA 4


Picture 32




Item No.	Description	Order-Code
	REPAIR KIT DISPLAY (EVITA 4) INCLUDING ALL PARTS AND THE LATEST DISPLAY VERSION	8411711
2	MULTIPLEXER NEC	8411686
3	LOOKUP TFT-TOSHIBA	8306795
6	SPACER	8300428
	LCD BACK LIGHTING (NEC), VERSION 2, without illustration	1841610

HINGED SUPPORT ARM

Picture 33




Item No.	Description	Order-Code
	HINGED SUPPORT ARM, W/O ILLUSTRATION, CURRENT VERSION W/O ITEM # 1 AND 2	DL-3020 MBL
1-7	HINGED ARM 240-DEGR.	8409609
1	SUPPORT FOR HINGED ARM	8409746
2	HOSE CLAMP FOR HINGED ARM	8409841
4	SET CLAMPING LEVER	8410383
5	BINDER, CPL.	8410075
6	STAR GRIP WITH DISC	8410382
7	RAIL CATCH	8411424


TROLLEY WITH 2 COLUMNS, SUPPORT**Picture 34**

Item No.	Description	Order-Code
1-2	SUPPORT	8412219
1	PLATE	M29110

TROLLEY WITH 2 COLUMNS, CARRIAGE

Picture 35





TROLLEY WITH 2 COLUMNS, BOTTOM SECTION**Picture 36**

Item No.	Description	Order-Code
1-12	TROLLEY	8411680
1	RECEIVER, CPL., SEE PICTURE 37 - 39	
2	CAP	8408853
3	SECTION TUBE	8408850
4	PROFILE	M29336
9a	SIDE PART, RIGHT, SEE PICTURE 40	
9b	SIDE PART, LEFT SEE PICTURE 41	
11	PIN	8408902

TROLLEY WITH 2 COLUMNS, RECEIVER


Picture 37




TROLLEY WITH 2 COLUMNS, SUPPORT**Picture 38**

TROLLEY WITH 2 COLUMNS, FLAP


Picture 39




Item No.	Description	Order-Code
1	DRAWER, CPL.	8409015
6	SHEET METAL	8408972
7	NUT	CH17646
8	SUPPORT	8408935
9	SLIDE	8408905
10	MAIN SPRING	M16452

TROLLEY WITH 2 COLUMNS, BOTTOM SECTION RIGHT**Picture 40**

Item No.	Description	Order-Code
1	FIXING GUIDE ROLL, ASSEMBLE CASTER WITH 1 STRAIN WASHER P/N 1341146	2M19140
5	GUIDE ROLL, ASSEMBLE CASTER WITH 1 STRAIN WASHER P/N 1341146	2M19130
	STRAIN WASHER DIN6796, w/o illustration	1341146

TROLLEY WITH 2 COLUMNS, BOTTOM SECTION LEFT**Picture 41**


Item No.	Description	Order-Code
1	FIXING GUIDE ROLL, ASSEMBLE CASTER WITH 1 STRAIN WASHER P/N 1341146	2M19140
5	GUIDE ROLL, ASSEMBLE CASTER WITH 1 STRAIN WASHER P/N 1341146	2M19130
	STRAIN WASHER DIN6796, w/o illustration	1341146


TROLLEY WITH 2 COLUMNS, FIXING SET FOR CUPBOARD**Picture 44**

Item No.	Description	Order-Code
1-2	CABINET ATTACHMENT SET EVITA	8409018


TROLLEY WITH 1 COLUMN

Picture 59




TROLLEY WITH 1 COLUMN**Picture 60**


Item No.	Description	Order-Code
1-17	MOBILE TROLLEY (HIGH)	8411950
2	COLUMN, HIGH, CPL.	8411952
4	COVER PLATE FOR UPPER PART	8411954
5	COVERING FOR BATTERIES	8411955
6	CASTOR WITH FIXING	8411945
7	CASTOR	8411946
8	HANDLE	8411957
9	HEADER PROFILE 25X5	8411958
10	SLIDE, CPL.	8411518
11	MAIN SPRING	M16452
12	SUPPORT	8408935
13	NUT	CH17646
14	CLOSING PLATE	8411959
15	SUPPORT, BACKSIDE	8411973
16	SHEET METAL	8411960

TROLLEY WITH 1 COLUMN, OPTION HUMIDIFIER RACK**Picture 61**

Item No.	Description	Order-Code
	HUMIDIFIER BRACKET EVITA MOBIL	8411956
5	SCREW M6X12 DIN 912	1263196

TROLLEY WITH 1 COLUMN, OPTION CYLINDER SUPPORT**Picture 62**


Item No.	Description	Order-Code
	CYLINDER SUPPORT, CPL., ORDER IN ADDITION HEIGHT EXTENSION P/N 8413962	8411970
1	CYLINDER SUPPORT	8411971
2	ADHESIVE CLOTH TAPE	8411972
	SCREW M6X12 DIN 912	1263196
	HEIGHT EXTENSION FOR TROLLEY, REQUIRED FOR CYLINDER SUPPORT	8413962

TROLLEY WITH 1 COLUMN, OPTION BATTERY CABLE**Picture 64**

Item No.	Description	Order-Code
1-4	KIT BATTERY CABLE EVITA MOBIL	8411822
	BATTERY 12V/17AH, ORDER 2 PIECES	1843303


FRONT PANEL EVITA 2 DURA

Picture 65



FRONT PANEL EVITA 2 DURA


Picture 66



Item No.	Description	Order-Code
1	CONTROL KNOB	M29954
2	SHAFT ENCODER, CPL.	8309118
3	KEYPAD	8411593
4	LOUDSPEAKER, CPL.	8306558
5	PCB 7-SEGMENT (EVITA 2 DURA)	8306871
6	PACKING RING	D04766
7	PCB FRONT PANEL (EVITA 2 DURA)	8306891
	COLOR TFT-DISPLAY, SEE PICTURE 68	

MODIFICATION KIT COLOR TFT-DISPLAY, EVITA 2 DURA

Picture 68



Item No.	Description	Part No.	Order-Code
1-4	KIT COLOUR PLUS		8411888
1	INVERTER SHARP		8411732
2	U DC/AC-CONVERTER 12V		1843648
3	CABLE CLIP		8712065
4	PACKING RING		D04766
	BACK LIGHT DISPLAY, SHARP, EVITA 2 DURA, WITHOUT ILLUSTRATION		1845357
	LCD BACK LIGHT, NEC, EVITA2 DURA, WITHOUT ILLUSTRATION		1844784

